
English	Instructions for use: Dental Implants	2
Español	Instrucciones de uso: Implantes Dentales	6
русский	Инструкция по применению: Зубные имплантаты	10
Deutsch	Gebrauchsanweisung: Dentalimplantate	14
Français	Instructions d'utilisation: Implants Dentaires	18
Italiano	Istruzioni per l'uso: Impianti Dentali	22
Português	Instruções de utilização: Implantes Dentários	26
Türkçe	Kullanım talimatları: Dental implantlar	30
简体中文	使用说明: 牙种植体	34
日本語	使用説明書: 歯科インプラント	37
한국어	적응증: 치과 시술용 임플란트	41
عربي	إرشادات الاستخدام: غرسات الأسنان	46
Polski	Instrukcja obsługi: Implanty dentystyczne	49
Українська	Інструкції з використання: зубні імплантати.	53
Čeština	Návod k použití: Zubní protetické komponenty	57

BioHorizons
2300 Riverchase Center
Birmingham AL, 35244 USA TOLL
FREE 888.246.8338
TEL 205.967.7880 FAX 205.870.0304
www.biohorizons.com

CE
2797

**ENGLISH**

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745
	Reference/ article number
	Lot/ batch number
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	EU Authorized Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Single sterile barrier system with protective packaging outside
	Single sterile barrier system
	Home

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Implants. Every product packaging label contains a description of the enclosed device.

DESCRIPTION

BioHorizons implants are manufactured from biocompatible titanium alloy and are single use only. Some BioHorizons implants come packaged with a 3inOne abutment pre-mounted. 3inOne abutments are for single patient use only. This abutment may be used as: (1) an implant mount; (2) a closed tray transfer when used with a Ball-top Screw and (3) a temporary or final abutment. "Peel-and-stick" labels supplied on the implant package contain important product information and should be applied to the patient's record in the event future reference is necessary.

INDICATIONS FOR USE

BioHorizons implants are intended for use in the mandible or maxilla for use as an artificial root structure for single tooth replacement or for fixed bridgework and dental retention. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion or (2) when splinted together for multiple tooth replacement or when stabilized with an overdenture supported by multiple implants. Tapered Short implants may be restored using delayed loading, or with a terminal or intermediate abutment for fixed or removable bridgework, and for overdentures. Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants may be used as an artificial root structure for single tooth replacement of mandibular central and lateral incisors and maxillary lateral incisors. The implants may be restored immediately (1) with a temporary prosthesis that is not in functional occlusion, (2) when splinted together as an artificial root structure for multiple tooth replacement of mandibular incisors, or (3) for denture stabilization using multiple implants in the anterior mandible and maxilla. The implants may be placed in immediate function when good primary stability has been achieved and with appropriate occlusal loading.

CONTRAINDICATIONS

BioHorizons dental implants should not be used in patients who have contraindicating systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants are contraindicated for single tooth replacement of maxillary central incisors.

DIRECTIONS FOR USE

Proper surgical procedures and techniques are the responsibility of the medical professional. Each dentist must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of postgraduate dental implant education and strict adherence to the instructions and procedures for use that accompany BioHorizons implant products. NOTE: A pre-operative 30-second rinse with a 0.12% Chlorhexidine Digluconate solution is recommended. (The influence of 0.12% Chlorhexidine Digluconate Rinses on the Incidence of Infectious Complications and Implant Success. Lambert, Paul M., et al, J Oral Maxillofacial Surgery 55:25-30, 1997, Suppl5). Tapered Internal Plus Implants - IMPORTANT! When using the TSK2021 (Tapered Internal Surgical Kit) to place the 3.8mm diameter implants, 3.0mm implant-level drivers are also required (TP3IDH, TP3IDR). Tapered Internal 3.0, Tapered Tissue Level 3.0, and Laser-Lok 3.0 Implants - IMPORTANT! (1) The implant procedure manual should be reviewed prior to attempting placement. (2) The osteotomy may be prepared with any BioHorizons 2.0mm and 2.5mm drill used in sequence. (3) Torque of 70Ncm or greater may compromise the mechanical integrity of the implant. In dense (D1-D2) bone further preparation of the site may be required to avoid applying excessive torque to the implant connection. (4) To further prepare the site, the crestal bone drill (TP3CBD) and the bone tap (122-900) may be used. (5) The clinician should verify that the proper instruments are available and at hand prior to attempting placement.

WARNINGS AND PRECAUTIONS

Appropriate training in proper implant surgery technique is strongly recommended prior to implant use. Improper technique can result in implant failure and/or loss of supporting bone. Appropriate x-ray films and/or CT scans should be utilized to determine (1) if adequate bone width and depth are available at the desired implant site and (2) the location of important anatomical landmarks, such as the mandibular canal, maxillary sinuses, and adjacent teeth. Clinician judgment should be used in determining the minimum post-implantation time before placing the implants in occlusal function. An adequate number of implants should be used to provide support and to distribute the load to the abutments. For short implants, clinicians should closely monitor patients for any of the following conditions: periimplant bone loss, changes to implant's response to percussion, or radiographic changes in bone to implant contact along the implant's length. If the implant shows mobility or greater than 50% bone loss, the implant should be evaluated for possible removal. If the clinicians choose a short implant, then clinicians should consider a two-stage surgical approach, splinting the short implant to an additional implant, and placement of the widest possible fixture. Allow longer periods for osseointegration avoid immediate loading. BioHorizons dental implants have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of BioHorizons dental implants in the MR environment is unknown. Scanning a patient who has this device may result in patient injury. Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. A certain percentage of implants may fail to achieve or to maintain osseointegration, as demonstrated by mobility, and should be removed. Potential causes of abutment fracture include, but not limited to: casting beyond a 30-degree angle, casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive cantilevering of pontics. Small diameter implants with either angled or straight abutments are intended for the anterior region of the mouth and not intended for the posterior region of the mouth due to possible failure of the implant. Tapered IM Implants (1) are limited to the molar region; (2) restricted to only one implant per molar extraction site; (3) are to be placed 3-5mm from adjacent implants or teeth, with the distance from the apical extent of the implant to adjacent implant or implant to tooth contact point should be 3mm to coronal to the alveolar crest; (4) should be at least 1.5mm from buccal or lingual alveolar surfaces to permit adequate integration; and, (5) are used when the surgeon has determined that the placement of a narrower diameter implant would increase the probability of failure due to poor primary stability, or increased surgical procedures leading to complications.

Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and competent authority of the EU Member State in which the clinician and/or patient is established.

COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) breakage of the implant required to be explanted and/or abutment required to be removed using clinician judgement; (3) loosening of the abutment screw and/or retaining screw; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) loosening of the implant requiring revision surgery; (9) perforation of the maxillary sinus; (10) perforation of the labial and lingual plates; and (11) bone loss possibly resulting in revision or removal.


HANDLING AND STERILIZATION

This product is supplied sterile and should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove from the package only after the correct size has been determined and the surgical site has been prepared. Always handle the product with powder-free gloves, and avoid contact with hard objects that may damage the surface. Implants are single use only, and re-use should not be attempted. Following this guidance eliminates the risk of cross-patient contamination from secondary use of this device. BioHorizons assumes no responsibility for attempted re-use or re-sterilization.

Products to be disposed of must be treated and decontaminated as dental surgery waste in compliance with the relevant regulations.

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