English	Instructions for use: Dental Prosthetics	2
Español	Instrucciones de uso: Prótesis Dentales	7
русский	Инструкция по применению: Зубные протезы	12
Deutsch	Gebrauchsanweisung: Dentale Prothetik	18
Français	Instructions d'utilisation: Prothèses Dentaires	23
Italiano	Istruzioni per l'uso: Protesi Dentali	28
Português	Instruções de utilização: Próteses Dentárias	33
Türkçe	Kullanım talimatları: Dental Protetik Ürünler	38
简体中文	使用说明: 牙科义齿	43
日本語	使用説明書: 歯科プロテーゼ	48
한국어	적응증: 치과 보철	54
عربي	إرشادات الاستخدام: البدليات السنية	60
Polski	Instrukcja obsługi: Protetyka Dentystyczna	65


BioHorizons  
2300 Riverchase Center  
Birmingham AL, 35244 USA TOLL  
FREE 888.246.8338  
TEL 205.967.7880 FAX 205.870.0304  
www.biohorizons.com

Quality First International OÜ  
Laki 30  
12915 Tallinn  
Estonia  
Tel +372-610-4196

**ENGLISH**

The symbol table below is for reference only. Refer to product packaging label for applicable symbols.

Symbol	Symbol Description
	Caution
	Electronic instructions for use
	Manufacturer
	BioHorizons products carrying the European Conformity (CE) mark fulfill the requirements of the Medical Device Directive 93/42/EEC as amended by Directive 2007/47/EC or the Medical Devices Regulation 2017/745
	Reference/ article number
	Lot/ batch number
	Do not re-use
	Do not re-sterilize
	Use-by-date
	Sterile by gamma irradiation
	Date of manufacture
	Caution: U.S. Federal law restricts these devices for sale, distribution and use by, or on the order of, a dentist or physician
	EU Authorized Representative
	Do not use if package is damaged
	Medical Device
	Non-sterile
	Home Position

This document supersedes all prior revisions. Original language is English.

This document applies to BioHorizons Dental Abutments and Prosthetic components. Every product packaging label contains a description of the enclosed device.

## **DESCRIPTION**

BioHorizons prosthetic components are intended for the restoration of BioHorizons dental implants within the specific indications of each implant system.

BioHorizons Multi-unit Abutments for CONELOG® and CONELOG OD Secure Abutments are intended for the restoration of CONELOG dental implants within the specific indications of each implant system.

The label on each prosthetic package contains important product information including whether the prosthetic is supplied sterile or non-sterile.

## **INDICATIONS FOR USE**

Tapered Internal and Internal Prosthetics are intended for the restoration of Tapered Internal and Internal dental implants and Zimmer® Dental ScrewVent® and Tapered ScrewVent® dental implants. Laser-Lok 3.0 Prosthetics are intended for the restoration of BioHorizons Laser-Lok 3.0 dental implants. Single-stage Prosthetics are intended for the restoration of BioHorizons Single-stage dental implants. External Prosthetics are intended for the restoration of BioHorizons External dental implants. Overdenture Prosthetics are intended for the restoration of BioHorizons Overdenture dental implants. One-piece Prosthetics are intended for the restoration of BioHorizons One-piece dental implants. Ceramic abutments are only intended for cement-retained cases in the anterior esthetic zone. BioHorizons Angled Multi-unit Abutments contain a unique LIGHT BLUE abutment screw specifically designed for use in those abutments. Multi-unit and Abutment for Screw prosthetics are intended for use as a support for multiple-implant supported prostheses. Ball Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible (3.0mm Ball Abutments being intended for use with 3.0mm dental implants). OD Secure Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla. Easy Ti Temp abutments are intended for long term temporary restorations (>30 days); hexed for single-unit and non-hexed for multiple-unit restorations. PEEK prosthetics are intended for short term temporary restorations (<30 days).

The BioHorizons Multi-unit Abutments for CONELOG are intended to function in the mandible or maxilla to support single and multiple-unit temporary or definitive restorations on CONELOG dental implants. CONELOG OD Secure Abutments are intended to retain a tissue-supported multiple-unit overdenture in the mandible and maxilla.

## **CONTRAINDICATIONS**

BioHorizons prosthetics should not be used in patients who have contraindicated systemic or uncontrolled local diseases such as blood dyscrasias, diabetes, hyperthyroidism, oral infections or malignancies, renal disease, uncontrolled hypertension, liver problems, leukemia, severe vascular heart disease, hepatitis, immunosuppressive disorder, pregnancy, collagen and bone diseases. Relative contraindications may include habits such as tobacco use, alcohol consumption, poor oral hygiene, bruxism, nail biting, pencil biting and improper tongue habits depending on severity. BioHorizons Multi-unit and Abutment for Screw prosthetics are contraindicated for single tooth replacements.

## **DIRECTIONS FOR USE**

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the Instructions

For Use that accompany BioHorizons products. Healing abutments should be hand tightened (10-15 Ncm). Abutment screws must be torqued to 30 Ncm to prevent screw loosening.

Platform	Zimmer Tapered ScrewVent Implants*	Zimmer ScrewVent Implants*
3.5mm Internal Connection	TSVBx TSV4Bx TSVT4Bx	TSVHx TSV4Hx
4.5mm Internal Connection	TSVWBx TSVTWBx TSVWHx	SVMBx SVBx SVMHx
5.7mm Internal Connection	TSV6Bx TSVT6Bx TSV6Hx	SVHx
		N/A

\*Where variable x = implant length.

Secure Custom Cast Abutments to the supplied polishing protector to protect the implant/abutment interface during laboratory procedures.

Single-stage Titanium Abutments for Laser Weld should be secured to the working model by torquing the accompanying screw to 30Ncm prior to laser welding. Secure Laser-Lok abutments onto the Laser-Lok Protective Sleeve to protect Laser-Lok during handling.

BioHorizons Multi-unit Abutment screws must be torqued to 30 Ncm.

The BioHorizons Multi-unit straight and angled abutment screws for CONELOG are packaged with the corresponding abutments and must be torqued to 20Ncm for intended function.

The BioHorizons Multi-unit Prosthetic Component screws must be torqued to 15Ncm.

BioHorizons Multi-unit Prosthetic Component design parameters are as follows:

- Minimum post height shall be equal to or greater than 4mm.
- Reduction of wall thickness is not permitted.
- Angle correction is not permitted. Angled Multi-unit abutments shall be used for angle correction.

BioHorizons Multi-unit Prosthetic Components used with the Multi-unit Abutments for CONELOG are intended to support multiple-unit temporary or definitive restorations except for the Titanium Coping. The Titanium Coping used with the Multi-unit Abutments for CONELOG is intended to support single and multiple-unit temporary or definitive restorations.

**WARNINGS AND PRECAUTIONS**

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Clinicians are responsible for understanding the appropriate technical use of BioHorizons prosthetic components. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at [www.biohorizons.com](http://www.biohorizons.com). Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU Member State in which the clinician and/or patient is established.

Dental implants can break in function for a number of reasons including overloading due to improper occlusion, metal fatigue, and over-tightening of the implant during insertion. Potential causes of abutment fracture include, but are not limited to: casting titanium above 2010°F (1099°C), inadequate implant support when attached to periodontically compromised teeth, non-passive fit of superstructure, overloading due to improper occlusion, incomplete seating of cemented abutments, and excessive

cantilevering of pontics. If any modifications are made to the implant/abutment interface, the abutment may not properly interface with the implant. The FDA considers the modifier of the implant/abutment interface a medical device company subject to FDA rules and regulations. Leaving bonding agent or particulate in the Laser-Lok surface may affect the connective tissue attachment. Use of the LIGHT BLUE abutment screws (packaged with Angled Multi-unit Abutments) with incompatible abutments may result in abutment and/or abutment screw failure.

Small diameter implants and angled abutments are not recommended for use in the posterior region of the mouth.

Prosthetics are single patient use only. To eliminate the risk of cross-patient contamination re-use should not be attempted. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

BioHorizons Multi-unit Abutments and Multi-unit Abutments for CONELOG require the addition of a coping or cylinder for single and multiple-unit restorations. The coping or cylinder gingival collar, angulation, and wall thickness shall not be modified. The coping or cylinder post height shall not be less than 4mm for single-unit restorations.

### **COMPLICATIONS AND ADVERSE EFFECTS**

The risks and complications with prosthetic components and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/ or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/ or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

### **HANDLING AND STERILIZATION**

Always handle the product with powder-free gloves and avoid contact with hard objects that may damage the surface. If the product is supplied sterile, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined and the surgical site has been prepared. Non-sterile product must be cleaned and sterilized. The following cleaning protocol must be used:

- 1) Prepare a detergent bath in a sterile container using a broad-spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® per the manufacturer's recommendations.
- 2) Brush the product to remove visible debris using a soft bristled brush moistened with the prepared detergent solution.
- 3) Thoroughly rinse product under running utility tap water.
- 4) Place product in the sterile container filled with the prepared detergent solution and sonicate for two (2) minutes minimum.
- 5) Thoroughly rinse product under running utility tap water.
- 6) Spray or wipe product with 70% IPA.
- 7) Blot product dry with clean lint free cloth.

For sterilization of non-sterile product (not including ceramic prosthetics), place product in an FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles:


Sterilization Cycles				
<b>Reference:</b>	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	ANSI/AAMI TIR12:2010	UK HTM 01-01 Part C:2016
<b>Type:</b>	Gravity Steam	Gravity Steam	Prevacuum Steam	Prevacuum Steam
<b>Exposure Time and Temperature:</b>	30 minutes at 121°C (250°F)	15 minutes at 132°C (270°F)	4 minutes at 132°C (270°F)	3 minutes at 134°C (273°F)
<b>Minimum Dry Time:</b>	30 minutes	30 minutes	20 minutes	20 minutes

For non-sterile ceramic prosthetics, dry heat sterilize at 170°C (338°F) for sixty (60) minutes.

It is recommended to include a 30-minute cool-down period before removing the product from the sterilization bag or wrap.

Zimmer®, Dental ScrewVent®, and Tapered ScrewVent® are registered trademarks of Zimmer Dental, Inc. Z-CAD® is a registered trademark of Metoxit AG. Conelog® is registered trademark of Camlog Biotechnologies GmbH.

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