## IMPORTANT MEDICAL INFORMATION DENTAL INSTRUMENTS

REV F APR 2016

Original Language is English

# BIOHORIZONS®

manufactured by:
**BioHorizons**
**2300 Riverchase Center Birmingham AL 35244 USA**
**TOLL FREE 888.246.8338**
**Tel. 205.967.7880 Fax 205.870.0304**
**www.biohorizons.com**

**EC REP**
**EU Authorised Representative**
**QUALITY FIRST INTERNATIONAL**
**Suites 317/318 Burford Business Centre 11 Burford Road, Stratford London E15 2ST, United Kingdom**
**Telephone +44-208-221-2361**
**Telefax +44-208-221-1912**

**CE**
**0086**

## ENGLISH

### DESCRIPTION

BioHorizons instruments are used for dental implant procedures such as site development, implant placement, and implant restorations within the specific indications of each implant system. The label on each instrument contains important product information including whether the instrument is supplied sterile or non-sterile and if the instrument is for single-use only. Instruments and kits supplied non-sterile must be cleaned and sterilized prior to first use and each use thereafter.

### INDICATIONS FOR USE

BioHorizons instruments and surgical kits are indicated for use in site development, placement, and restoration of BioHorizons implants and associated components. The BioHorizons Single Use Drills are indicated for use with the BioHorizons Tapered Internal and Laser-Lok 3.0 implants only. No attempt should be made to re-use or re-sterilize the drills. BioHorizons assumes no responsibility for attempted re-use or re-sterilization between patients.

### CONTRAINDICATIONS

BioHorizons instruments should not be used with patients that have allergies to the specific materials used, including stainless steel and titanium alloy.

### DIRECTIONS FOR USE

Proper surgical procedures and restorative techniques are the responsibility of the medical professional. Each clinician must evaluate the appropriateness of the procedure used based on personal medical training and experience as applied to the patient case at hand. BioHorizons strongly recommends completion of dental implant courses and strict adherence to the instructions pertaining to BioHorizons products.

**Internal Abutment-level Driver, Hex-chuck Handpiece:** Use only with compatible W&H Hexagon Chucking System handpieces. **Drills and Taps:** Should be replaced when wear is noticed, such as a decrease in cutting efficiency or when signs of discoloration appear. Cutting instruments should be replaced after approximately 12 to 20 osteotomy cycles, depending on bone density. BioHorizons recommends use of a drill usage chart to track drill use and to ensure drills are replaced as directed. **Single Use Drills:** Insert the drill into a handpiece and drills protective tubing off prior to use. **CGS Drills and Guides:** Care must be taken to ensure the CGS Drills follow the angle of the guide handle inserts or sleeves to prevent drill binding. Refer to the drill report for proper drill selection and usage. **Site Dilators:** Push the dilator into the bone while rotating it. In firmer bone, it may be necessary to use a mallet to tap the dilator to proper depth. **Laser-Lok Protective Sleeves:** Secure Titanium Base Abutments onto the sleeve side labeled Hybrid or Custom Titanium Abutments onto the side labeled Full Ti to protect Laser-Lok region during bonding and polishing.

### WARNINGS AND PRECAUTIONS

Clinician judgment, as related to individual patient presentations, must always supersede recommendations in any BioHorizons Instructions for Use (IFU). Refer to OSHA standard 29CFR1910.1030 prior to cleaning and sterilization. Additional technical information is available upon request from BioHorizons, or may be viewed and/or downloaded at www.biohorizons.com. Contact BioHorizons Customer Care or your local representative with any questions you have regarding a specific IFU.

**Prosthetic Instruments** - Only use a torque wrench for final installation of prosthetic components requiring a specified installation torque. Finger-tightening of prosthetic components can result in insufficient torque and eventual loosening of the component. **Over-tightening prosthetic components could break the component or spin the mating implant.**

### COMPLICATIONS AND ADVERSE EFFECTS

The risks and complications with instruments, prosthetics, and implants include, but are not limited to: (1) allergic reaction(s) to implant and/or abutment material; (2) implant and/or abutment breakage; (3) abutment screw and/or retaining screw loosening; (4) infection requiring revision of the dental implant; (5) nerve damage that could cause permanent weakness, numbness, or pain; (6) histologic responses possibly involving macrophages and/or fibroblasts; (7) formation of fat emboli; (8) implant loosening requiring revision surgery; (9) maxillary sinus perforation; (10) labial or lingual plate perforation; and (11) bone loss possibly resulting in revision or removal.

### HANDLING AND STERILIZATION

Always handle the product with power-free gloves and avoid contact with hard objects that may damage the surface. If the **product is supplied sterile**, it should be considered sterile unless the package has been opened or damaged. Using accepted sterile technique, remove product from the package only after the correct size has been determined. If the **product is supplied non-sterile or when reprocessing**: Remove and discard any shipping material before initial sterilization. The device must be thoroughly cleaned. Remove any visible debris from the instruments and surgical

kit using a soft bristle brush and a broad spectrum cleaning or disinfecting agent such as Hu-Friedy's Enzymax® or equivalent. Rinse thoroughly. Place the instruments in a beaker of the same solution and sonicate for 10 minutes. Rinse thoroughly. Rinse instruments with isopropyl alcohol to remove any soap residue and minerals. Blow instruments with a lint free towel and air dry completely. **For AS123 Hand Unit:** blow pressurized air through one of the four drain holes to purge. **CGS Quick Connect Handle:** blow pressurized air through one end to purge. **For both instruments:** use a medical instrument lubricant such as STERIS® Hinge Free® after each use.) Refer to the labeling of the cleaning/lubricating agent used for instructions for use. Instruments that are unable to be cleaned, discolored, do not properly interface with mating components, and/or do not articulate as designed should be disposed of. Disassemble the surgical kit and wash the empty tray using the same solution. Rinse with water and dry thoroughly. Return the instruments to the appropriate locations in the tray. Place product in a FDA cleared sterilization bag or wrap and run through one of the following qualified sterilization cycles: 1. Prevacuum Steam: 132°C (270°F) for five (5) minutes minimum. 2. Gravity Steam: 132°C (270°F) for thirty (30) minutes minimum. 3. Gravity Steam: 121°C (250°F) for sixty (60) minutes minimum. 4. Prevacuum Steam: 134°C (273°F) for three (3) minutes minimum. Dry 20 to 50 minutes. Do not exceed 134°C (273°F) during sterilization or drying.

**Attention!** Proper cleaning may lead to inadequate sterilization. Failure to completely dry instruments during autoclaving may leave moisture and cause discoloration and oxidation. The use of hydrogen peroxide or other oxidizing agents will damage the surface of the instruments. Periodic testing, cleaning, and calibration of the autoclave equipment is recommended to ensure the unit remains in proper working order.

CAUTION: FEDERAL LAW (US) RESTRICTS THIS DEVICE TO THE SALE, DISTRIBUTION, AND USE BY OR ON THE ORDER OF A LICENSED MEDICAL PROFESSIONAL.

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