

Includes the instrumentation required to place: Tapered Pro, Tapered Plus, Tapered Tissue Level, Tapered 3.0 and Tapered Internal.

2.0mm HD Drill TSD2020HD

2.5mm HD Drill TSD2025HD

2.8mm HD Drill TSD2028HD

3.2mm HD Drill TSD2032HD

3.7mm HD Drill TSD2037HD

4.1mm HD Drill TSD2041HD

4.7mm HD Drill TSD2047HD

5.4mm HD Drill TSD2054HD

3.0mm HD Depth Gauge TDG2030HD

3.4mm HD Depth Gauge TDG2034HD

3.8mm HD Depth Gauge TDG2038HD

4.6mm HD Depth Gauge TDG2046HD

5.8mm HD Depth Gauge TDG2058HD

3.0mm HD Crestal Bone Drill TSC2030HD

3.4mm HD Crestal Bone Drill TSC2034HD

3.8mm HD Crestal Bone Drill TSC2038HD

4.6mm HD Crestal Bone Drill TSC2046HD

5.8mm HD Crestal Bone Drill TSC2058HD

3.0mm HD Bone Tap TST2030HD

3.4mm HD Bone Tap TST2034HD

3.8mm HD Bone Tap TST2038HD

4.6mm HD Bone Tap TST2046HD

5.8mm HD Bone Tap TST2058HD

Drill Extender 122-100

4mm Square Drive Extender 300-205

Hand Wrench 300-400

.050" (1.25mm) Hex Driver 135-351

Drill Sequence (Color = Implant Diameter)

5.8, 4.6, 3.8, 3.4, 3.0

2.0, 2.5, 3.2, 3.7, 4.1, 4.7, 5.4

Drivers (Color = Prosthetic Platform)

0, 7.5, 9, 10.5, 12, 15, 18

Handpiece, Implant, Abutment, Ratchet

BIOHORIZONS® tapered internal HD

Straight Parallel Pins
(2 per kit)
144-100

20° Angled Parallel Pins
(2 per kit)
144-200

3.0 Implant-level Driver, Handpiece, Regular
TP3IDHR

3.5/4.5mm HD Implant-level Driver, Handpiece
TYGIDH

5.7mm HD Implant-level Driver, Handpiece
TBIDH

Internal Abutment-level Driver, Handpiece
PHA

3.0 Implant-level Driver, Ratchet, Regular
TP3IDRR

3.5/4.5mm HD Implant-level Driver, Ratchet
TYGIDR

5.7mm HD Implant-level Driver, Ratchet
TBIDR

Internal Abutment-level Driver, Ratchet
PRA

Ratchet
130-000

Spare slots can be used to hold preferred surgeon instrumentation such as stop drills and extended shank drills.

The implant staging area is designed to temporarily secure implant vials during surgery only. Do not attempt to sterilize implant vials in the surgical tray.

customer care: 888.246.8338 • shop online at store.biohorizons.com



Tapered Internal HD
surgical kit (TSK4000)



L02025



REV D FEB 2020